


 <p><b>62</b> W6: 41,2 LC: 8,4 WC: 7,7</p>	 <p><b>1N</b> W6: 43,6 LC: 10,0 WC: 8,0</p>
 <p><b>3M</b> W6: 40,3 LC: 9,0 WC: 7,3</p>	 <p><b>2N</b> W6: 42,1 LC: 10,6 WC: 7,9</p>
 <p><b>2D</b> W6: 42,7 LC: 9,6 WC: 7,9</p>	 <p><b>3N</b> W6: 44,4 LC: 10,3 WC: 8,3</p>
 <p><b>3D</b> W6: 44,5 LC: 9,5 WC: 8,4</p>	 <p><b>2P</b> W6: 46,0 LC: 11,2 WC: 8,3</p>
 <p><b>90</b> W6: 44,1 LC: 9,4 WC: 8,0</p>	 <p><b>3P</b> W6: 47,6 LC: 10,7 WC: 8,8</p>
 <p><b>15</b> W6: 43,6 LC: 10,3 WC: 8,0</p>	 <p><b>21</b> W6: 45,3 LC: 10,9 WC: 8,5</p>
 <p><b>A25</b> W6: 44,7 LC: 9,6 WC: 8,5</p>	 <p><b>1H</b> W6: 48,2 LC: 11,9 WC: 8,8</p>
 <p><b>A26</b> W6: 47,6 LC: 10,2 WC: 9,2</p>	 <p><b>20</b> W6: 50,1 LC: 12,5 WC: 9,4</p>
 <p><b>4H</b> W6: 47,4 LC: 10,1 WC: 8,8</p>	 <p><b>23</b> W6: 49,5 LC: 13,3 WC: 9,4</p>
 <p><b>136</b> W6: 49,1 LC: 10,4 WC: 9,1</p>	<p>*Humedeciendo el talón de los dientes con monómero, se fortalece la unión con la base de la dentadura.                      * إن ترطيب خلفية الأسنان بمادة أحادية الخلايا، يقوي من التلاحم مع قاعدة الأسنان.                      *Base of moistened with a monomer strengthens the bond to denture base.</p>



	<b>4</b> W6: 33,5 LC: 8,4 WC: 5,0
	<b>5</b> W6: 32,2 LC: 8,9 WC: 4,6
	<b>9</b> W6: 34,7 LC: 8,3 WC: 5,0
	<b>2DD</b> W6: 31,9 LC: 8,8 WC: 4,9
	<b>3PP</b> W6: 37,1 LC: 9,8 WC: 5,7
	<b>1NN</b> W6: 33,1 LC: 9,0 WC: 5,1
	<b>3NN</b> W6: 34,9 LC: 9,8 WC: 5,2
	<b>13</b> W6: 36,8 LC: 9,7 WC: 5,3
	<b>226</b> W6: 35,0 LC: 10,6 WC: 5,3
	<b>3R</b> W6: 36,6 LC: 11,2 WC: 5,3

	 <b>50</b> 30 32,7 8,9 9,3
	 <b>52</b> 31,4 35,1 9,2 10,0
	 <b>54</b> 33,9 38,1 9,7 10,4
	 <b>29</b> 28,9 31,6 7,5 8,8
	 <b>31</b> 31,6 34,2 8,1 9,2
	 <b>32</b> 31,4 34,6 8,7 9,3





ARTICULACIONES ACONSEJABLES  
RECOMENDED ARTICULATIONS  
الفكوك المنصوح بها  
ARTICULATIONS CONSEILLÉES

ANTERIORES ANTERIORES الفك الأمامي ANRIEURES			POSTERIORES POSTERIORES الفك الخلفي POSTERIEURES		
SUP. الفك العلوي SUP.	INFERIORES INFERIOR الفك السفلي INFERIEURES				
62	4		50	29	29MZ
3M	5	4	50	31	29MZ
2D	2DD	4	50	31	29MZ
3D	3PP	3NN	50	31	29MZ
90	9		52	32	33M
15	1NN		52	33	33M
A25	3PP		52	32	33M
A26	3PP		54	33	33M
4H	3PP		54	33	33M
136	13		54	33	33M
1N	3NN	1NN	54	32	33L
2N	1NN		52	32	33M
3N	3NN		52	32	33M
2P	226		54	34	33L
3P	3PP		54	33	33M
21	13		54	33	33L
1H	226		54	34	33L
20	3R		54	34	
23	3R		54	34	

في حالات العضات القصيرة يستحسن استعمال الموديلات الخاصة.  
For short cases use specially moulds.

NORMAL		SHORT BITE	
33°	50-52-54-34	20°	29LL
20°	29-31-33-32		
0°	29MZ-33M-33L	10°	F32

LAS CIFRAS SIGNIFICAN: MEDIDAS (mm)  
THE FIGURES ARE MEANING: MEASURES (mm)  
الأرقام تعني: مقاييس  
LES NOMBRES SIGNIFIANT: MESURES (mm)

